

# MEDICAL INFRARED THERMOMETER

RAPID INDUCTION • ACCURATE, NONCONTACT MEASUREMENT



a key  
measure



forehead  
measurement



memory  
storage



hot tip



replaceable  
battery



LED  
display



a second  
display



automatic  
shutdown



c/f  
switch



intelligent  
accurate

This medical grade device can be used as part of your required COVID Health Assessment Program for employees and customers. Keep your employees safe, maintain business continuity and build consumer confidence with proper health screening, alongside social distancing, sanitizing and handwashing programs.

\*Color of product may vary.  
Distributed by DayMark Safety Systems.

THE DIFFERENCE IS NIGHT AND DAY.™

800.847.0101

daymark.com

#3517B-DIS  
121724

**DayMark**<sup>®</sup>  
Safety Systems


# MEDICAL INFRARED THERMOMETER

RAPID INDUCTION • ACCURATE, NONCONTACT MEASUREMENT




This highly sensitive probe measures temperature in as little as one second. Suitable for quick measurement of multiple team members or taking an individual's temperature multiple times.

## Unit Description



## Display Information



## Technical Specifications

**Measurement Mode:** Forehead

**Operating Conditions:** 60.8° F to 95° F (16° C to 35° C)  
with relative humidity of 15% to 80%

**Measuring Range of Body Temperature:** 89.6° F to 108° F  
(32° C - 42.2° C)

**Accuracy:** ± .54° F (± .3° C)

**Measurement Distance:** .5" to 1.95" (1.5 cm to 5 cm) from  
the measuring point

**Resolution:** 32.18° F (.1° C)

**Memory:** 32 Measurement Recalls

**Power Supply:** (2) AAA (DC 3V) batteries

**Dimension:** 5.57" x 3" x 1.5" (144mm x 78mm x 37mm)

**Weight:** 3.1 oz. (88 g) excludes battery

**Auto Power Off:** 60 seconds

**FDA, CE & CFDA Certified**  
**RoHS Report Available**

\*Color of product may vary.  
Distributed by DayMark Safety Systems.

THE DIFFERENCE IS NIGHT AND DAY.™

800.847.0101

daymark.com

**DayMark®**  
Safety Systems